



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## BHLHB3 Antibody [#C11661]

**Catalog Number:** C11661

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9C0J9

**Other Names:** bHLHB3, Differentially expressed in chondrocytes protein 2, hDEC2, Enhancer-of-split and hairy-related protein 1, SHARP-1

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human BHLHB3.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** BHLHB3 antibody detects endogenous levels of total BHLHB3 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:40000

**References:**

Fujimoto K., Biochem. Biophys. Res. Commun. 280:164-171(2001).

Garriga-Canut M., J. Biol. Chem. 276:14821-14828(2001).

