



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ILKAP Antibody [#C11670]

Catalog Number: C11670

Concentration: 1mg/ml

Swiss-Prot No.: Q9H0C8

Other Names: EC 3.1.3.16, ILKAP, Integrin-linked kinase-associated serine/threonine phosphatase 2C

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ILKAP.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ILKAP antibody detects endogenous levels of total ILKAP protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:20000

References:

Leung-Hagesteijn C., EMBO J. 20:2160-2170(2001).


Wiemann S., Genome Res. 11:422-435(2001).

Hillier L.W., Nature 434:724-731(2005).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue using ILKAP antibody.



Peptide - +

Western blot analysis of extracts from COS-7 cells, using ILKAP antibody.