



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ZNF541 Antibody [#C11671]

Catalog Number: C11671

Concentration: 1mg/ml

Swiss-Prot No.: Q9H0D2

Other Names: Zinc finger protein 541, ZNF541, ZN541

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ZNF541.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ZNF541 antibody detects endogenous levels of total ZNF541 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Wiemann S., Genome Res. 11:422-435(2001).

Yu L.-R., J. Proteome Res. 6:4150-4162(2007).

Sjoebloom T., Science 314:268-274(2006).

