



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ECRG4 Antibody [#C11679]

Catalog Number: C11679

Concentration: 1mg/ml

Swiss-Prot No.: Q9H1Z8

Other Names: Esophageal cancer-related gene 4 protein
[Precursor], C2orf40, ECRG4, UNQ761/PRO1508

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ECRG4.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity: ECRG4 antibody detects endogenous levels of total ECRG4 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)

Applications: IHC: 1:50~1:100 IF: 1:100~1:500
ELISA: 1:40000


References:

Bi M., Submitted (DEC-2000) to the EMBL/GenBank/DDBJ databases.
Clark H.F., Genome Res. 13:2265-2270(2003).
Hillier L.W., Nature 434:724-731(2005).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human testis tissue, using ECRG4 antibody.



Peptide - +

Immunofluorescence analysis of A549 cells, using ECRG4 antibody.