



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

ADNP Antibody [#C11685]

Catalog Number: C11685

Concentration: 1mg/ml

Swiss-Prot No.: Q9H2P0

Other Names: Activity-dependent neuroprotective protein, Activity-dependent neuroprotector, ADNP, ADNP1, KIAA0784

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

ADNP.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ADNP antibody detects endogenous levels of total ADNP protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Zamostiano R., J. Biol. Chem. 276:708-714(2001).

Deloukas P., Nature 414:865-871(2001)

Nagase T., DNA Res. 5:277-286(1998).

