



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

HTR5B Antibody [#G016]

Catalog Number: G016

Concentration: 1mg/ml

Swiss-Prot No.: P35365

Other Names: 5-hydroxytryptamine receptor 5B; 5-HT-5B;
Serotonin receptor 5B; MR22; 5ht5b

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human HTR5B.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: HTR5B antibody detects endogenous levels of total HTR5B protein.


Reactivity: Mouse, Rat

Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:40000


References:

Erlander M.G., Proc. Natl. Acad. Sci. U.S.A. 90:3452-3456(1993).

Wisden W., FEBS Lett. 333:25-31(1993).



Western blot analysis of extracts from COLO cells, using HTR5B antibody.



Immunofluorescence analysis of HUVEC cells, using HTR5B antibody.