



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

ADORA3 Antibody [#G021]

Catalog Number: G021

Concentration: 1mg/ml

Swiss-Prot No.: P33765

Other Names: A3AR; AA3R; adenosine A3 receptor

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ADORA3.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ADORA3 antibody detects endogenous levels of total ADORA3 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:5000

References:


Sajjadi F.G., *Biochim. Biophys. Acta* 1179:105-107(1993).

Salvatore C.A., *Proc. Natl. Acad. Sci. U.S.A.* 90:10365-10369(1993).

Gregory S.G., *Nature* 441:315-321(2006).



Western blot analysis of extracts from HepG2 cells, using ADORA3 antibody.



Immunofluorescence analysis of COS-7 cells, using ADORA3 antibody.