



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

ADRA2A Antibody [#G022]

Catalog Number: G022

Concentration: 1mg/ml

Swiss-Prot No.: P08913

Other Names: Alpha-2A adrenergic receptor; Alpha-2A adrenoceptor; Alpha-2A adrenoreceptor; Alpha-2AAR; Alpha-2 adrenergic receptor subtype C10; ADRA2A; ADRA2R, ADRAR

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

ADRA2A.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ADRA2A antibody detects endogenous levels of total ADRA2A protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:20000

References:

Kobilka B.K., Science 238:650-656(1987).

Fraser C.M., J. Biol. Chem. 264:11754-11761(1989).

Guyer C.A., J. Biol. Chem. 265:17307-17317(1990).

