



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** +86 519 8805 0026

## ADRA1B Antibody [#G024]

**Catalog Number:** G024

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P35368

**Other Names:** Alpha-1B adrenergic receptor; Alpha 1B-adrenoceptor; Alpha 1B-adrenoreceptor; ADRA1B

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human ADRA1B.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ADRA1B antibody detects endogenous levels of total ADRA1B protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%);  
Mouse (Identities = 93%, Positives = 93%);  
Rat (Identities = 93%, Positives = 93%)

**Applications:** IF: 1:100~1:500 ELISA: 1:20000

**References:**

Ramarao C.S., J. Biol. Chem. 267:21936-21945(1992).

Forray C., Mol. Pharmacol. 45:703-708(1994).

Schwinn D.A., J. Pharmacol. Exp. Ther. 272:134-142(1995).

