



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

AGTRL1 Antibody [#G031]

Catalog Number: G031

Concentration: 1mg/ml

Swiss-Prot No.: P35414

Other Names: Apelin receptor; G-protein coupled receptor APJ; Angiotensin receptor-like 1; HG11; AGTRL1; APJ

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human AGTRL1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AGTRL1 antibody detects endogenous levels of total AGTRL1 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000

References:


O'Dowd B.F., Gene 136:355-360(1993).

Eggerickx D., Submitted (JUN-1995) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team; Genome Res. 14:2121-2127(2004).



Western blot analysis of extracts from 293 cells, using AGTRL1 antibody.



Immunofluorescence analysis of MCF-7 cells, using AGTRL1 antibody.