



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: +86 519 8805 0026

GPR132 Antibody [#G130]

Catalog Number: G130

Concentration: 1mg/ml

Swiss-Prot No.: Q9UNW8

Other Names: Probable G-protein coupled receptor 132; G2 accumulation protein; GPR132; G2A

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GPR132.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GPR132 antibody detects endogenous levels of total GPR132 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Weng Z., Proc. Natl. Acad. Sci. U.S.A. 95:12334-12339(1998).

Kabarowski J.H.S., Science 293:702-705(2001).

Witte O.N., Science 307:206-206(2005).

