



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

CHRM4 Antibody [#G208]

Catalog Number: G208

Concentration: 1mg/ml

Swiss-Prot No.: P08173

Other Names: Muscarinic acetylcholine receptor M4;
CHRM4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CHRM4.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CHRM4 antibody detects endogenous levels of total CHRM4 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%);
Rat (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:10000

References:

Bonner T.I., Science 237:527-532(1987).

Bonner T.I., Neuron 1:403-410(1988).

Peralta E.G., EMBO J. 6:3923-3929(1987).

