



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPR171 Antibody [#G143]

Catalog Number: G143

Concentration: 1mg/ml

Swiss-Prot No.: O14626

Other Names: Probable G-protein coupled receptor 171;
G-protein coupled receptor H963; GPR171; H963

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human GPR171.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: GPR171 antibody detects endogenous levels of total GPR171 protein.


Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:1000


References:

Jacobs K.A., Gene 198:289-296(1997).

Joensuu T., Am. J. Hum. Genet. 69:673-684(2001).



Western blot analysis of extracts from K562 cells, using GPR171 antibody.



Immunofluorescence analysis of MCF-7 cells, using GPR171 antibody.