



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

GPR173 Antibody [#G144]

Catalog Number: G144

Concentration: 1mg/ml

Swiss-Prot No.: Q9NS66

Other Names: Probable G-protein coupled receptor 173;
Super conserved receptor expressed in brain 3; GPR173;
SREB3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from internal of human
GPR173.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: GPR173 antibody detects endogenous levels
of total GPR173 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Matsumoto M., Biochem. Biophys. Res. Commun.
272:576-582(2000).

The MGC Project Team; Genome Res.
14:2121-2127(2004).

