



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

ADORA1 Antibody [#G201]

Catalog Number: G201

Concentration: 1mg/ml

Swiss-Prot No.: P30542

Other Names: Adenosine receptor A1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human ADORA1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ADORA1 antibody detects endogenous levels of total ADORA1 protein.

Reactivity: Human, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Libert F., Biochem. Biophys. Res. Commun. 187:919-926(1992).

Townsend-Nicholson A., Brain Res. Mol. Brain Res. 16:365-370(1992).

Ren H., J. Biol. Chem. 269:3104-3110(1994).

