



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

ADORA3 Antibody [#G204]

Catalog Number: G204

Concentration: 1mg/ml

Swiss-Prot No.: P33765

Other Names: Adenosine A3 receptor;ADORA3;
UNQ1931/PRO4406

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ADORA3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Specificity: ADORA3 antibody detects endogenous levels of total ADORA3 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 92%, Positives = 92%);
Rat (Identities = 92%, Positives = 92%)

Applications: IF: 1:100~1:500 ELISA: 1:10000

References:

Sajjadi F.G., *Biochim. Biophys. Acta* 1179:105-107(1993).
Salvatore C.A., *Proc. Natl. Acad. Sci. U.S.A.* 90:10365-10369(1993).
Meyerehof W., Submitted (DEC-1993) to the EMBL/GenBank/DDBJ databases.

