



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

CHRM1 Antibody [#G205]

Catalog Number: G205

Concentration: 1mg/ml

Swiss-Prot No.: P11229

Other Names: ACM1; muscarinic acetylcholine receptor M1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CHRM1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CHRM1 antibody detects endogenous levels of total CHRM1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Allard W.J., Nucleic Acids Res. 15:10604-10604(1987).

Chapman C.G., Nucleic Acids Res. 18:2191-2191(1990).

Peralta E.G., EMBO J. 6:3923-3929(1987).

