



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TAF4 Antibody [#C10018]

Catalog Number: C10018

Concentration: 1mg/ml

Swiss-Prot No.: O00268

Other Names: Transcription initiation factor TFIID subunit 4, TBP-associated factor 4, Transcription initiation factor TFIID 135 kDa subunit, TAF(II)135, TAFII-135, TAFII135, TAFII-130, TAFII130, RNA polymerase II TBP-associated factor subunit C, TAF4, TAF2C, TAF2C1, TAF4A, TAFII130, TAFII135

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from internal of human TAF4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TAF4 antibody detects endogenous levels of total TAF4 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Mengus G., Genes Dev. 11:1381-1395(1997).

Deloukas P., Nature 414:865-871(2001).

Tanese N., Proc. Natl. Acad. Sci. U.S.A. 93:13611-13616(1996).



For Research Use Only