



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TR11B Antibody [#C10019]

Catalog Number: C10019

Concentration: 1mg/ml

Swiss-Prot No.: O00300

Other Names: Tumor necrosis factor receptor superfamily member 11B [Precursor], Osteoprotegerin, Osteoclastogenesis inhibitory factor, TNFRSF11B, OCIF, OPG

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from Internal of human

TR11B.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TR11B antibody detects endogenous levels of total TR11B protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Simonet W.S., Cell 89:309-319(1997).

Yasuda H., Endocrinology 139:1329-1337(1998).

Morinaga T., Eur. J. Biochem. 254:685-691(1998).

