



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

LHR2A Antibody [#C10028]

Catalog Number: C10028

Concentration: 1mg/ml

Swiss-Prot No.: O00534

Other Names: Loss of heterozygosity 11 chromosomal region 2 gene A protein, Breast cancer suppressor candidate 1, BCSC-1, LOH11CR2A, BCSC1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from Internal of human LHR2A.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: LHR2A antibody detects endogenous levels of total LHR2A protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Martin E.S., Proc. Natl. Acad. Sci. U.S.A. 100:11517-11522(2003).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Monaco C., Genomics 46:217-222(1997).

