



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

DLGP1 Antibody [#C10034]

Catalog Number: C10034

Concentration: 1mg/ml

Swiss-Prot No.: O14490

Other Names: DAP-1, DLGAP1, DLGP1, Disks

large-associated protein 1, Guanylate kinase-associated protein, PSD- 95/SAP90-binding protein 1, SAP90/PSD-95-associated protein 1, SAPAP1, hGKAP

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

DLGP1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: DLGP1 antibody detects endogenous levels of total DLGP1 protein.


Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:5000


References:

Satoh K., Genes Cells 2:415-424(1997).

Kim E., J. Cell Biol. 136:669-678(1997).



Western blot analysis of extracts from HepG2 cells, using DLGP1 antibody.



Immunofluorescence analysis of HeLa cells, using DLGP1 antibody.