



Web: [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: (86) 519 8805 0026

## TP53I11 Antibody [#C10041]

**Catalog Number:** C10041

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O14683

**Other Names:** Tumor protein p53-inducible protein 11, p53-induced gene 11 protein, TP53I11, PIG11

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

TP53I11.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** TP53I11 antibody detects endogenous levels of total TP53I11 protein.


**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 IHC: 1:50~1:100

ELISA: 1:20000


**References:**

Polyak K., Nature 389:300-306(1997).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human brain tissue using TP53I11 antibody.



TP53I11 --

Peptide - +

Western blot analysis of extracts from HUVEC cells, using TP53I11 antibody.