



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

BACH1 Antibody [#C10055]

Catalog Number: C10055

Concentration: 1mg/ml

Swiss-Prot No.: O14867

Other Names: BAC1, BTB and CNC homolog 1, HA2303, Transcription regulator protein BACH1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human BACH1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: BACH1 antibody detects endogenous levels of total BACH1 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:20000

References:


Blouin J.-L., Hum. Genet. 102:282-288(1998).

Ohira M., Genomics 47:300-306(1998).

Hattori M., Nature 405:311-319(2000).



Western blot analysis of extracts from HeLa cells, HepG2 cells and 293 cells, using BACH1 antibody.



Immunofluorescence analysis of HeLa cells, using BACH1 antibody.