



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

KLF11 Antibody [#C10058]

Catalog Number: C10058

Concentration: 1mg/ml

Swiss-Prot No.: O14901

Other Names: KLF11, Krueppel-like factor 11, TGFB-inducible early growth response protein 2, TIEG-2, transforming growth factor-beta-inducible early growth response protein 2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from Internal of human

KLF11.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: KLF11 antibody detects endogenous levels of total KLF11 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Cook T., J. Biol. Chem. 273:25929-25936(1998).

Asano H., Mol. Cell. Biol. 19:3571-3579(1999).

Neve B., Proc. Natl. Acad. Sci. U.S.A. 102:4807-4812(2005).

