



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

GAK Antibody [#C10064]

Catalog Number: C10064

Concentration: 1mg/ml

Swiss-Prot No.: O14976

Other Names: Cyclin G-associated kinase, EC 2.7.11.1,
kinase GAK

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from Internal of human GAK.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: GAK antibody detects endogenous levels of
total GAK protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Kimura S.H., Genomics 44:179-187(1997).

The MGC Project Team; Genome Res.
14:2121-2127(2004).

Greener T., J. Biol. Chem. 275:1365-1370(2000).

