



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TPD54 Antibody [#C10100]

Catalog Number: C10100

Concentration: 1mg/ml

Swiss-Prot No.: O43399

Other Names: hD54; TPD54; tumor protein D52-like 2; tumor protein D54

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human TPD54.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TPD54 antibody detects endogenous levels of total TPD54 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Nourse C.R., Biochim. Biophys. Acta 1443:155-168(1998).

Deloukas P., Nature 414:865-871(2001).

Wilson S.H.D., Genomics 76:81-88(2001).

