



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AN32C Antibody [#C10101]

Catalog Number: C10101

Concentration: 1mg/ml

Swiss-Prot No.: O43423

Other Names: Acidic leucine-rich nuclear phosphoprotein 32 family member C, Phosphoprotein 32-related protein 1, Tumorigenic protein pp32r1, ANP32C, PP32R1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

AN32C.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AN32C antibody detects endogenous levels of total AN32C protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Kadkol S.S., Nat. Med. 5:275-279(1999).

Kochevar G.J., Hum. Mutat. 23:546-551(2004).

Matilla A., Cerebellum 4:7-18(2005).

