



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TBX1 Antibody [#C10102]

Catalog Number: C10102

Concentration: 1mg/ml

Swiss-Prot No.: O43435

Other Names: CAFS, CTHM, DGCR, DGS, DORV, T-box 1, T-box 1 transcription factor C

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human TBX1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TBX1 antibody detects endogenous levels of total TBX1 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:20000

References:


Chieffo C., Genomics 43:267-277(1997).

Yagi H., Lancet 362:1366-1373(2003).

Sjoeblom T., Science 314:268-274(2006).



Western blot analysis of extracts from COLO205 cells, using TBX1 antibody.



Immunofluorescence analysis of A549 cells, using TBX1 antibody.