



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AKAP10 Antibody [#C10110]

Catalog Number: C10110

Concentration: 1mg/ml

Swiss-Prot No.: O43572

Other Names: A kinase anchor protein 10, mitochondrial [Precursor], Protein kinase A-anchoring protein 10, PRKA10, Dual specificity A kinase-anchoring protein 2, D-AKAP-2, AKAP10

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from Internal of human AKAP10.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AKAP10 antibody detects endogenous levels of total AKAP10 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Wang L., Proc. Natl. Acad. Sci. U.S.A. 98:3220-3225(2001).

Chatterjee T.K., Submitted (DEC-1997) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team, Genome Res. 14:2121-2127(2004).

