



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

HAIR Antibody [#C10111]

Catalog Number: C10111

Concentration: 1mg/ml

Swiss-Prot No.: O43593

Other Names: Protein hairless, HR, HAIR

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human HAIR.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: HAIR antibody detects endogenous levels of total HAIR protein.

Reactivity: Human (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 IF: 1:100~1:500
ELISA: 1:10000

References:

Ahmad W., Science 279:720-724(1998).


Cichon S., Hum. Mol. Genet. 7:1671-1679(1998).

Ahmad W., Am. J. Hum. Genet. 63:984-991(1998).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using HAIR antibody.



Peptide - +

Immunofluorescence analysis of A549 cells, using HAIR antibody.