



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## EFNA2 Antibody [#C10124]

**Catalog Number:** C10124

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O43921

**Other Names:** Ephrin-A2 [Precursor], EPH-related receptor tyrosine kinase ligand 6, LERK-6, HEK7-ligand, HEK7-L, EFNA2, EPLG6, LERK6

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from Internal of human EFNA2.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** EFNA2 antibody detects endogenous levels of total EFNA2 protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:5000

**References:**

Cerretti D.P., Genomics 47:131-135(1998).

Aasheim H.C., Biochem. Biophys. Res. Commun. 252:378-382(1998).

Grimwood J., Nature 428:529-535(2004).

