



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ZNF387 Antibody [#C10130]

Catalog Number: C10130

Concentration: 1mg/ml

Swiss-Prot No.: O60284

Other Names: Suppression of tumorigenicity protein 18, Zinc finger protein 387, ST18, KIAA0535, ZNF387

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human ZNF387.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ZNF387 antibody detects endogenous levels of total ZNF387 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Nagase T., DNA Res. 5:31-39(1998).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Jandrig B., Oncogene 23:9295-9302(2004).

