



Web: www.anobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

HBP1 Antibody [#C10135]

Catalog Number: C10135

Concentration: 1mg/ml

Swiss-Prot No.: O60381

Other Names: HMG box-containing protein 1, HMG box-containing protein 1, High mobility group box transcription factor 1, HBP1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from Internal of human

HBP1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: HBP1 antibody detects endogenous levels of total HBP1 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Lemercier C., Mol. Cell. Biol. 20:6627-6637(2000).

Hillier L.W., Nature 424:157-164(2003).

Ota T., Nat. Genet. 36:40-45(2004).



For Research Use Only