



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

hnRNP Q Antibody [#C10139]

Catalog Number: C10139

Concentration: 1mg/ml

Swiss-Prot No.: O60506

Other Names: Heterogeneous nuclear ribonucleoprotein Q, hnRNP Q, hnRNP-Q, Synaptotagmin-binding, cytoplasmic RNA-interacting protein, Glycine- and tyrosine-rich RNA-binding protein, GRY-RBP, NS1-associated protein 1, SYNCRIP, HNRPQ, NSAP1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from Internal of human

hnRNP Q.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: hnRNP Q antibody detects endogenous levels of total hnRNP Q protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:20000

References:


Harris C.E., J. Virol. 73:72-80(1999).

Du G., Chin. Sci. Bull. 45:343-349(2000).

Mourelatos Z., EMBO J. 20:5443-5452(2001).



Western blot analysis of extracts from HUVEC cells, Jurkat cells and 293 cells, using hnRNP Q antibody.



Immunofluorescence analysis of HepG2 cells, using hnRNP Q antibody.

For Research Use Only