



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

SDCG1 Antibody [#C10140]

Catalog Number: C10140

Concentration: 1mg/ml

Swiss-Prot No.: O60524

Other Names: NY-CO-1, SDCG1, serologically defined colon cancer antigen 1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human SDCG1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: SDCG1 antibody detects endogenous levels of total SDCG1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Ota T., Nat. Genet. 36:40-45(2004).

Scanlan M.J., Int. J. Cancer 76:652-658(1998).

Olsen J.V., Cell 127:635-648(2006).

