



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TRIM24 Antibody [#C10073]

Catalog Number: C10073

Concentration: 1mg/ml

Swiss-Prot No.: O15164

Other Names: TIF1, TIF1A, TRIM24, Transcription intermediary factor 1-alpha, Tripartite motif protein 24, kinase TIF1-alpha

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human TRIM24.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TRIM24 antibody detects endogenous levels of total TRIM24 protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:20000

References:

Thenot S., J. Biol. Chem. 272:12062-12068(1997).

Venturini L., Oncogene 18:1209-1217(1999).


Hellal-Levy C., Mol. Endocrinol. 14:1210-1221(2000).



Peptide -

+

Immunohistochemistry analysis of paraffin-embedded human brain tissue using TRIM24 antibody.



Peptide -

+

Western blot analysis of extracts from COLO cells, using TRIM24 antibody.