



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

p73 Antibody [#C10080]

Catalog Number: C10080

Concentration: 1mg/ml

Swiss-Prot No.: O15350

Other Names: p53-like transcription factor; p53-related protein; P73 alpha protein; p73-alpha; TP73; TRP73; Tumor protein p73

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human p73.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: p73 antibody detects endogenous levels of total p73 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:40000

References:

Kaghad M., Cell 90:809-819(1997).


Mai M., Genomics 51:359-363(1998).

De Laurenzi V., J. Exp. Med. 188:1763-1768(1998).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using p73 antibody.



Peptide - +

Western blot analysis of extracts from Jurkat cells, using p73 antibody.

For Research Use Only