



Web: [www.anbobio.com](http://www.anbobio.com)

[www.dijibio.com](http://www.dijibio.com)

Order: [order@dijibio.com](mailto:order@dijibio.com)

Tel: (86) 519 8805 0026

## M3K13 Antibody [#C10092]

**Catalog Number:** C10092

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O43283

**Other Names:** Mitogen-activated protein kinase kinase kinase 13, EC 2.7.11.25, Mixed lineage kinase, MLK, Leucine zipper-bearing kinase, MAP3K13, LZK

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

M3K13.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** M3K13 antibody detects endogenous levels of total M3K13 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Sakuma H., J. Biol. Chem. 272:28622-28629(1997).

Schultz S.J., Cell Growth Differ. 4:821-830(1993).

Ikeda A., FEBS Lett. 488:190-195(2001).

