



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

LMX1B Antibody [#C10143]

Catalog Number: C10143

Concentration: 1mg/ml

Swiss-Prot No.: O60663

Other Names: LIM homeobox transcription factor 1 beta, LIM/homeobox protein LMX1B, LIM/homeobox protein 1.2, LMX-1.2, LMX1B

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

LMX1B.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: LMX1B antibody detects endogenous levels of total LMX1B protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Dreyer S.D., Nat. Genet. 19:47-50(1998).

Vollrath D., Hum. Mol. Genet. 7:1091-1098(1998).

McIntosh I., Am. J. Hum. Genet. 63:1651-1658(1998).

