



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

DNA Polymerase ζ Antibody [#C10144]

Catalog Number: C10144

Concentration: 1mg/ml

Swiss-Prot No.: O60673

Other Names: DNA polymerase zeta catalytic subunit, hREV3, REV3L, POLZ, REV3

Storage/Stability: Store at -20°C /1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human DNA Polymerase ζ .

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: DNA Polymerase ζ antibody detects endogenous levels of total DNA Polymerase ζ protein.

Reactivity: Human, Mouse


Applications: IHC: 1:50~1:100 ELISA: 1:20000

References:

Gibbs P.E.M., Proc. Natl. Acad. Sci. U.S.A. 95:6876-6880(1998).

Lin W., Mutat. Res. 433:89-98(1999).

Morelli C., Cytogenet. Cell Genet. 83:18-20(1998).



Peptide - +
Immunohistochemistry analysis of paraffin-embedded human brain tissue using DNA Polymerase ζ antibody.