



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

BRS3 Antibody [#G217]

Catalog Number: G217

Concentration: 1mg/ml

Swiss-Prot No.: P32247

Other Names: Bombesin receptor subtype-3; BRS-3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human BRS3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: BRS3 antibody detects endogenous levels of total BRS3 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Fathi Z., J. Biol. Chem. 268:5979-5984(1993).

Gorbulev V., FEBS Lett. 340:260-264(1994).

Ross M.T., Nature 434:325-337(2005).

