



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

CCBP2 Antibody [#G221]

Catalog Number: G221

Concentration: 1mg/ml

Swiss-Prot No.: O00590

Other Names: Chemokine-binding protein 2; Chemokine-binding protein D6; C-C chemokine receptor D6; Chemokine receptor CCR-9; Chemokine receptor CCR-10; CCBP2; CCR10, CMKBR9

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

CCBP2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CCBP2 antibody detects endogenous levels of total CCBP2 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:20000

References:

Bonini J.A., DNA Cell Biol. 16:1249-1256(1997).

Nibbs R.J.B., J. Biol. Chem. 272:32078-32083(1997).

Kopatz S.A., Submitted (MAR-2003) to the EMBL/GenBank/DDBJ databases.

