



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## PFKFB2 Antibody [#C10151]

**Catalog Number:** C10151

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O60825

**Other Names:** 6-phosphofructo-2-kinase,  
6-phosphofructo-2-kinase/fructose-2,6-biphosphatase 2,  
6PF-2-K/Fru-2,6-P2ase heart-type isozyme, F262,  
PFK-2/FBPase-2, PFK-2/FBPase-2, heart type,  
fructose-2,6-bisphosphatas

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered  
saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl,  
0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against

synthesized peptide derived from N-terminal of human  
PFKFB2.

**Purification:** The antibody was affinity-purified from  
rabbit antiserum by affinity-chromatography using  
epitope-specific immunogen.

**Specificity:** PFKFB2 detects endogenous levels of total  
PFKFB2 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Heine-Suner D., Eur. J. Biochem. 254:103-110(1998).

Soejima H., Genomics 74:115-120(2001).

Gregory S.G., Nature 441:315-321(2006).

