



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

CCKAR Antibody [#G222]

Catalog Number: G222

Concentration: 1mg/ml

Swiss-Prot No.: P32238

Other Names: CCK-A receptor; CCK-AR; CCKR; CCKRA;
cholecystokinin type A receptor

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CCKAR.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CCKAR antibody detects endogenous levels of total CCKAR protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Ulrich C.D., Biochem. Biophys. Res. Commun. 193:204-211(1993).

Wank S.A., Biochem. Biophys. Res. Commun. 194:811-818(1993).

Miller L.J., Gastroenterology 109:1375-1380(1995).

