



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## Dyskerin Antibody [#C10152]

**Catalog Number:** C10152

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O60832

**Other Names:** CBF5 homolog, DKC1, EC 5.4.99.-,  
Nucleolar protein NAP57

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from Internal of human Dyskerin.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** Dyskerin antibody detects endogenous levels of total Dyskerin protein.

**Reactivity:** Human, Mouse, Rat


**Applications:** WB: 1:500~1:3000 IF: 1:100~1:500  
ELISA: 1:10000

**References:**


Heiss N.S., Nat. Genet. 19:32-38(1998).

Knight S.W., Am. J. Hum. Genet. 65:50-58(1999).

Hassock S., Genomics 55:21-27(1999).



Western blot analysis of extracts from HeLa cells, using Dyskerin antibody.



Immunofluorescence analysis of HeLa cells, using Dyskerin antibody.