



Web: [www.anbobio.com](http://www.anbobio.com)

[www.dijibio.com](http://www.dijibio.com)

Order: [order@dijibio.com](mailto:order@dijibio.com)

Tel: (86) 519 8805 0026

## TBX10 Antibody [#C10164]

**Catalog Number:** C10164

**Concentration:** 1mg/ml

**Swiss-Prot No.:** O75333

**Other Names:** T-box transcription factor TBX10, T-box protein 10, TBX10, TBX7

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human TBX10.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** TBX10 antibody detects endogenous levels of total TBX10 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:5000

**References:**

Bush J.O., Gene Expr. Patterns 3:533-538(2003).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Law D.J., Mamm. Genome 9:397-399(1998).

