



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

EMR4 Antibody [#G246]

Catalog Number: G246

Concentration: 1mg/ml

Swiss-Prot No.: Q86SQ3

Other Names: egf-like module containing, mucin-like, hormone receptor-like 4; emr4; gpr127; pgr16

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human EMR4.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: EMR4 antibody detects endogenous levels of total EMR4 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Fredriksson R., Biochem. Biophys. Res. Commun. 301:725-734(2003).

Hamann J., Eur. J. Immunol. 33:1365-1371(2003).

Vassilatis D.K., Proc. Natl. Acad. Sci. U.S.A. 100:4903-4908(2003).

