



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

FFAR2 Antibody [#G249]

Catalog Number: G249

Concentration: 1mg/ml

Swiss-Prot No.: O15552

Other Names: Free fatty acid receptor 2; G-protein coupled receptor 43; FFAR2; GPR43

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human FFAR2.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FFAR2 antibody detects endogenous levels of total FFAR2 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:1000

References:

Sawzdargo M., Biochem. Biophys. Res. Commun. 239:543-547(1997).

Grimwood J., Nature 428:529-535(2004).

Kaighin V.A., Submitted (DEC-2007) to the EMBL/GenBank/DDBJ databases.

