



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TACC1 Antibody [#C10172]

Catalog Number: C10172

Concentration: 1mg/ml

Swiss-Prot No.: O75410

Other Names: gastric cancer antigen Ga55; KIAA1103; taxin 1; transforming acidic coiled-coil-containing protein 1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human TACC1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TACC1 antibody detects endogenous levels of total TACC1 protein.

Reactivity: Human, Mouse


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000

References:


Still I.H., *Oncogene* 18:4032-4038(1999).

Still I.H., *BMC Evol. Biol.* 4:16-16(2004).

Kikuno R., *DNA Res.* 6:197-205(1999).



Western blot analysis of extracts from K562 cells, using TACC1 antibody.



Immunofluorescence analysis of MCF-7 cells, using TACC1 antibody.