



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MED24 Antibody [#C10176]

Catalog Number: C10176

Concentration: 1mg/ml

Swiss-Prot No.: O75448

Other Names: 100 kDa thyroid hormone receptor associated protein, ARC100, Activator-recruited cofactor 100 kDa component, DRIP100, KIAA0130, THRAP4, TR100, Thyroid hormone receptor-associated protein complex 100 kDa component, Vitamin D3 receptor-interacting protein complex component DRIP100

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

MED24.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MED24 antibody detects endogenous levels of total MED24 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:20000

References:


Yuan C.-X., Proc. Natl. Acad. Sci. U.S.A. 95:7939-7944(1998).

Rachez C., Genes Dev. 12:1787-1800(1998).

Nagase T., DNA Res. 2:167-174(1995).



Western blot analysis of extracts from 293 cells, using MED24 antibody.



Immunofluorescence analysis of HUVEC cells, using MED24 antibody.